



LOCKING  
IMPLANTS

 **SAHA**<sup>®</sup>  
Orthopedic Implants

Manufacturer of  
**Orthopedic** Implants &  
Surgical **Instruments**

## SMALL LOCKING COMPRESSION PLATE

Product Code: **B0.155**

<b>Holes</b>	3 Hole to 16 Hole
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## SMALL T OBLIQUE LCP – LEFT & RIGHT

Product Code: **B0.151**

<b>Holes</b>	2 Hole to 10 Hole
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium

## SMALL 'T' BUTTRESS LCP

Product Code: **B0.150**

<b>Holes</b>	2 Hole to 10 Hole
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## RECONSTRUCTION LCP (Dynamic Hole) 3.5MM

Product Code: B0.9912

<b>Holes</b>	3 Hole to 16 Hole
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## SUPERIOR CLAVICLE LCP – LEFT & RIGHT

Product Code: B0.9928

<b>Holes</b>	4H, 6H, 8H
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium

## CLAVICLE PLATE - II LCP LEFT & RIGHT

Product Code: B0.01.9928

<b>Holes</b>	6 Hole, 8 Hole, 10 Hole
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium



## ONE THIRD TUBULAR LCP ROUND HOLE 3.5MM

Product Code: **B0.01.121**

<b>Holes</b>	3 Hole to 16 Hole
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## ONE THIRD TUBULAR LCP DC HOLE 3.5MM

Product Code: **B0I.03.121**

<b>Holes</b>	4Hole to 14 Hole
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium

## LATERAL END CLAVICLE LCP (2.7MM & 3.5MM) L & R

Product Code: **B0.501**

<b>Holes</b>	3 Hole To 8 Hole
<b>Locking Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



**PROXIMAL HUMERUS LCP PHILOS  
(LONG DYNAMIC HOLE)**

Product Code: **B0I.156**

<b>Holes</b>	2H,3H4H,5H,6H,7H,8H,9H,10H,11H,12H
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Locking Screw 4.0mm</b>	20mm to 70mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



**PROXIMAL HUMERUS LCP PHILOS  
(SMALL DYNAMIC HOLE)**

Product Code: **B0.02.156**

<b>Holes</b>	2 Hole to 13 Hole
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium

**PROXIMAL HUMERUS LCP**

Product Code: **B0.9911**

<b>Holes</b>	3H,4H,5H,6H,7H,8H,9H,10H
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Locking Screw 4.0mm</b>	20mm to 70mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## CLAVICLE HOOK LCP 3.5MM – LEFT & RIGHT

Product Code: **B0.270**

<b>Hook Depths</b>	12mm,15mm,18mm
<b>Holes</b>	4H,5H,6H,7H
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## CLOVERLEAF LCP

Product Code: **B0.154**

<b>Holes</b>	4H,5H,6H,7H,8H,9H,10H,11H,12H
<b>Locking Screw 3.5mm</b>	10mm to 50mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium

## MEDIAL DISTAL TIBIA LCP – LEFT & RIGHT

Product Code: **B0.9909**

<b>Holes</b>	5H,7H,9H,11H,13H,15H
<b>Locking Screw 3.5mm</b>	10mm to 50mm (Diff. 2.0mm)
<b>Locking Screw 4.0mm</b>	20mm to 70mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## DISTAL TIBIA ANATOMICAL LCP – LEFT & RIGHT

Product Code: **B0.2001**

<b>Holes</b>	4H,6H,8H,10H,12H,14H
<b>Locking Screw 3.5mm</b>	10mm to 50mm (Diff. 2.0mm)
<b>Locking Screw 4.0mm</b>	20mm to 70mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## CALCANEAL LCP 3.5MM - LEFT & RIGHT

Product Code: **B0.9940**

<b>Holes</b>	65MM, 70MM, 75MM
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium

## DISTAL TIBIA LCP WITH TAB 3.5MM - LEFT & RIGHT

Product Code: **B0.267**

<b>Holes</b>	5H,7H,9H,11H,13H
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Locking Screw 4.0mm</b>	20mm to 70mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## DISTAL FIBULA LCP

Product Code: **B0.269**

<b>Holes</b>	3 Hole To 10 Hole
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## DISTAL FIBULA PLATE LCP (DLX) – LEFT & RIGHT

Product Code: **B0.01.269**

<b>Holes</b>	3 Hole to 10 Hole
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium

## DISTAL FIBULA LCP 2.7MM X 3.5MM

Product Code: **B0.02.269**

<b>Holes</b>	4 Hole To 10 Hole
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Cortical Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium



## RECONSTRUCTION LCP 5MM (DYNAMIC HOLE)

Product Code: **B0.9913**

<b>Holes</b>	4H,5H,6H,7H,8H,9H,10H,11H,16H
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## NARROW LOCKING COMPRESSION PLATE

Product Code: **B0.166**

<b>Holes</b>	6H,7H,8H,9H,10H,11H,12H,16H
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium

## BROAD LOCKING COMPRESSION PLATE

Product Code: **B0.167**

<b>Holes</b>	6H,7H,8H,9H,10H,11H,12H,16H
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## 'L' BUTTRESS LCP 5MM - LEFT & RIGHT

Product Code: **B0.148**

<b>Holes</b>	4H,5H,6H,7H,8H,9H,10H,11H,12H
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## T LOCKING COMPRESSION PLATE 5MM

Product Code: **B0.162**

<b>Holes</b>	4H,5H,6H,7H,8H,9H,10H,11H,12H
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium

## 'T' BUTTRESS LCP 5MM

Product Code: **B0I.262**

<b>Holes</b>	4H,5H,6H,7H,8H,9H,10H,11H,12H
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## LATERAL TIBIAL LCP 6.5MM - LEFT & RIGHT

Product Code: **B0.161**

<b>Holes</b>	4H,5H,6H,7H,8H,9H,10H,11H,12H
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Locking Screw 6.5mm</b>	35mm to 90mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## LATERAL TIBIAL LCP 5.0MM - LEFT & RIGHT

Product Code: **B0.159**

<b>Holes</b>	5H,7H,9H,11H,13H,15H
<b>Cancellous Locking Screw 5.0mm</b>	30mm to 80mm (Diff. 5.0mm)
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## MEDIAL PROXIMAL TIBIA LCP 3.5MM - LEFT & RIGHT

Product Code: **B0.268**

<b>Holes</b>	4H,6H,8H,10H,12H,14H,16H
<b>Locking Screw 3.5mm</b>	10mm to 50mm (Diff. 2.0mm)
<b>Locking Screw 4.0mm</b>	20mm to 85mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## POSTEROMEDIAL PROXIMAL TIBIA LCP 3.5MM

Product Code: **B0.273**

<b>Holes</b>	3H,4H,5H,6H,7H,8H,9H,10H,11H,12H
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Locking Screw 4.0mm</b>	20mm to 85mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## DISTAL HUMERUS EXTRA-ARTICULAR LCP – L & R

Product Code: **B0.9941**

<b>Holes</b>	4H,6H,8H,10H,12H,14H
<b>Locking Screw 3.5mm</b>	10mm to 50mm (Diff. 2.0mm)
<b>Locking Screw 4.0mm</b>	20mm to 70mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium

## MEDIAL DISTAL HUMERUS LCP (2.7MM & 3.5MM) L & R

Product Code: **B0.9908**

<b>Holes</b>	3H,5H,7H,9H,11H,13H
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Locking Screw 2.7mm</b>	12mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## MEDIAL DISTAL FEMUR LCP 5MM – LEFT & RIGHT

Product Code: **B0.231**

<b>Holes</b>	4H,5H,6H
<b>Locking Screw 5.0mm</b>	20mm to 50mm (Diff. 2.0mm)
	55mm to 95mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## DISTAL HUMERUS LCP (2.7MM & 3.5MM) – LEFT & RIGHT

Product Code: **B0.9925**

<b>Holes</b>	3H,5H,7H,9H,11H,13H
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Locking Screw 2.7mm</b>	20mm to 40mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium

## DISTAL HUMERUS WITH LATERAL SUPPORT LCP (2.7MM & 3.5MM) – LEFT & RIGHT

Product Code: **B0.9907**

<b>Holes</b>	3H,5H,7H,9H,11H,13H
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Locking Screw 2.7mm</b>	12mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## DISTAL FEMUR LCP 6.5MM - LEFT & RIGHT

Product Code: **B0.163**

<b>Holes</b>	5H,7H,9H,11H,13H,15H
<b>Cancellous Locking Screw 6.5mm</b>	30mm to 95mm (Diff. 5.0mm)
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## DISTAL FEMUR LCP 5MM - LEFT & RIGHT

Product Code: **B0.164**

<b>Holes</b>	5H,7H,9H,11H,13H,15H
<b>Cancellous Locking Screw 5.0mm</b>	30mm to 80mm (Diff. 5.0mm)
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium

## DHS LCP (LONG & SHORT BARREL)

Product Code: **B0.165 & B0.152**

<b>Holes</b>	3H,4H,5H,6H,7H,8H,9H,10H
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Lag Screw</b>	50mm to 115mm (Diff. 5.0mm)
<b>Angel</b>	130° , 135°
<b>Material</b>	Stainless Steel, Titanium



## 95° D.C.S LCP

Product Code: **B0.9937**

<b>Holes</b>	6H,7H,8H,9H,10H,11H,12,13H,14H,15H,16H
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Cancellous Locking Screw 6.5mm</b>	30mm to 95mm (Diff. 5.0mm)
<b>Lag Screw</b>	50mm to 115mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## PROXIMAL FEMUR LCP COBRA HEAD (8MM SCREW)

Product Code: **B0.03.9938**

<b>Holes</b>	5H,7H,9H,11H,13H,15H
<b>Cannulated Screw 8.0mm</b>	55mm to 120mm (Diff. 5.0mm)
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Cancellous Locking Screw 6.5mm</b>	30mm to 95mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium

## PROXIMAL FEMUR LCP COBRA HEAD 6.5MM

Product Code: **B0.02.9938**

<b>Holes</b>	3H, 5H, 7H, 9H, 11H, 13Hole
<b>Locking Screw 6.5mm</b>	30mm to 95mm (Diff. 5.0mm)
<b>Cancellous Screw 6.5mm</b>	30mm to 95mm (Diff. 5.0mm)
<b>Material</b>	Stainless steel, Titanium



## PROXIMAL FEMUR LCP

Product Code: **B0.9938**

<b>Holes</b>	5H,7H,9H,11H,13H,15H
<b>Locking Screw 5.0mm</b>	10mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Cancellous Screw 6.5mm</b>	30mm to 100mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## PROXIMAL FEMUR LCP 6.5MM

Product Code: **B0.01.9938**

<b>Holes</b>	3H, 5H, 7H, 9H, 11H, 13Hole
<b>Locking Screw 6.5mm</b>	30mm to 95mm (Diff. 5.0mm)
<b>Material</b>	Stainless steel, Titanium

## PROXIMAL LATERAL TIBIA LCP 3.5MM - LEFT & RIGHT

Product Code: **B0.9986**

<b>Holes</b>	4H,6H,8H,10H,12H,14H
<b>Locking Screw 3.5mm</b>	12mm to 40mm (Diff. 2.0mm)
<b>Locking Screw 4.0mm</b>	20mm to 85mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## PROXIMAL TIBIA LOCKING PLATE NEW 3.5MM LEFT & RIGHT

Product Code: B0.01.9986

<b>Holes</b>	4H, 6H, 8H, 10H, 12H, 14Hole
<b>Locking Screw 3.5mm</b>	10mm to 60mm (Diff. 2.0mm)
<b>Cortical Screw 3.5mm</b>	10mm to 60mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium



## TOMIFIX MEDIAL HIGH TIBIAL LCP

Product Code: B0I.252

<b>Holes</b>	3H,4H
<b>Locking Screw 5.0mm</b>	18mm to 50mm (Diff. 2.0mm)
	55mm to 90mm (Diff. 5.0mm)
<b>Cancellous Locking Screw 5.0mm</b>	35mm to 100mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## DISTAL RADIUS LCP - STRAIGHT

Product Code: B0I.193

<b>Holes</b>	4H,5H,6H
<b>Locking Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.7mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## 'L' DISTAL RADIUS LCP (HEAD 2 HOLE) – LEFT & RIGHT

Product Code: **B0I.191**

<b>Holes</b>	3H,4H
<b>Locking Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.7mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## 'L' DISTAL RADIUS LCP (HEAD 3HOLE) – LEFT & RIGHT

Product Code: **B0I.190**

<b>Holes</b>	3H,4H
<b>Locking Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.7mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## 'L' DISTAL RADIUS OBLIQUE LCP (HEAD 3HOLE) – L & R

Product Code: **B0I.194**

<b>Holes</b>	3H,4H
<b>Locking Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.7mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## 'T' DISTAL RADIUS LCP (HEAD 3HOLE)

Product Code: **B0I.192**

<b>Holes</b>	3H,4H
<b>Locking Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.7mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## DISTAL RADIUS LCP (HEAD 5HOLE) – LEFT & RIGHT

Product Code: **B0I.197**

<b>Holes</b>	3H,4H,5H,6H,7H
<b>Locking Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.7mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## DISTAL RADIUS BUTTRESS LCP (HEAD 5HOLE) – L & R

Product Code: **B0I.198**

<b>Holes</b>	3Hole
<b>Locking Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.7mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



**TWO-COLUMA DISTAL RADIUS PLATE – LEFT & RIGHT**

Product Code: **B0.9946**

<b>Holes</b>	2H,3H,4H,5H,6H,7H,8H
<b>Locking Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.7mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



**VOLAR DISTAL RADIUS (DLX) LCP PLATE STANDARD L & R**

Product Code: **B0.272**

<b>Holes</b>	2H,3H,4H,5H,6H,7H,8H
<b>Locking Screw 3.5mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Locking Screw 2.7mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.7mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



**VOLAR DISTAL RADIUS (DLX) LCP WIDE – LEFT & RIGHT**

Product Code: **B0.01.272**

<b>Holes</b>	2H,3H,4H,5H,6H,7H,8H
<b>Locking Screw 3.5mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Locking Screw 2.7mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.7mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## RADIAL HEAD LCP RIM PLATE 2.4MM

Product Code: **B0.182**

<b>Holes</b>	2H,3H,4H,5H,6H,7H
<b>Locking Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## RADIAL HEAD LCP 2.4MM – LEFT & RIGHT

Product Code: **B0.183**

<b>Holes</b>	2H,3H,4H,5H,6H,7H
<b>Locking Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Cortex Screw 2.4mm</b>	06mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium

## OLECRANON LCP 3.5MM - LEFT & RIGHT

Product Code: **B0I.509**

<b>Holes</b>	4H,6H,8H,10H
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## ANTEROLATERAL DISTAL TIBIA LCP 3.5MM L & R

Product Code: **B0.511**

<b>Holes</b>	5H,7H,9H,11H,13H,15H
<b>Locking Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Locking Screw 4.0mm</b>	20mm to 70mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## PROXIMAL FEMUR LOCKING PLATE (PEDIATRIC)


Product Code: **B0.9939**

<b>Holes</b>	3Hole, 4Hole, 5Hole
<b>Degree</b>	110, 130, 140, 150
<b>Locking Screw 3.5mm</b>	10mm to 60mm (Diff. 2.0mm)
<b>Cortical Screw 3.5mm</b>	10mm to 60mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium

## DUAL LOCK-PROXIMAL FEMUR LOCKING PLATE 6.0MM LEFT & RIGHT

Product Code: **B0.9943**

<b>Holes</b>	4 Hole, 6 Hole, 8 Hole, 10 Hole, 12 Hole, 14 Hole
<b>Cannulated Locking Screw 6.0mm</b>	50mm to 120mm (Diff. 5.0mm)
<b>Material</b>	Stainless steel, Titanium



## ANKLE AND FOOT BONE (FIRST MTP FUSION PLATE 5<sup>0</sup>) LOCKING PLATE - LEFT & RIGHT

Product Code: **B0.531**

<b>Holes</b>	Small, Medium
<b>Locking Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Cortical Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Titanium



## ANKLE AND FOOT BONE (FIRST TMT FUSION) LOCKING PLATE

Product Code: **B0.532**

<b>Holes</b>	Standard, Long
<b>Locking Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Cortical Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Titanium

## ANKLE AND FOOT BONE (SECOND TMT FUSION) LOCKING PLATE

Product Code: **B0I.533**

<b>Locking Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Cortical Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Titanium



## ANKLE AND FOOT BONE (CLOVERLEAF FUSION) LOCKING PLATE

Product Code: B0.534

<b>Holes</b>	Short, Standard
<b>Locking Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Cortical Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Titanium



## ANKLE AND FOOT BONE (T FUSION) LOCKING PLATE

Product Code: B0.535

<b>Holes</b>	Short, Standard
<b>Locking Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Cortical Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Titanium

## ANKLE AND FOOT BONE (L FUSION) LOCKING PLATE - LEFT & RIGHT

Product Code: B0.536

<b>Holes</b>	Short, Standard
<b>Locking Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Cortical Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Titanium



## ANKLE AND FOOT BONE (X) LOCKING PLATE

Product Code: **B0.537**

<b>Holes</b>	Extra Small, Small, Medium, Large
<b>Locking Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Cortical Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Titanium



## ANKLE AND FOOT BONE (STRAIGHT FUSION) LOCKING PLATE

Product Code: **B0.538**

<b>Holes</b>	4 Hole
<b>Locking Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Cortical Screw 2.7mm</b>	10mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Titanium

## ANKLE LOCKING PLATE 3.5MM

Product Code: **B0.514**

<b>Holes</b>	4 Hole
<b>Locking Screw 3.5mm</b>	10mm to 60mm (Diff. 2.0mm)
<b>Cortical Screw 3.5mm</b>	10mm to 60mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium



### MICRO RECON PLATE 2.0MM

Product Code: B0.521

<b>Holes</b>	4 Hole, 6 Hole, 8 Hole, 10 Hole, 12 Hole
<b>Cortical Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Material</b>	Titanium



### MICRO CONDYLAR LOCKING PLATE

Product Code: B0.522

<b>Holes</b>	3 Hole, 5 Hole, 7 Hole
<b>Locking Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Cortical Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Material</b>	Titanium

### MICRO Y LOCKING PLATE 2.0MM

Product Code: B0.523

<b>Holes</b>	3 Hole, 5 Hole, 7 Hole
<b>Locking Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Cortical Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Material</b>	Titanium



### MICRO T LOCKING PLATE 2.0MM

Product Code: B0.524

<b>Holes</b>	3 Hole, 5 Hole, 7 Hole
<b>Locking Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Cortical Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Material</b>	Titanium



### MICRO ROTATION CORRECTION LOCKING PLATE LCP 2.0MM

Product Code: B0.525

<b>Holes</b>	4 Hole, 5 Hole
<b>Locking Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Cortical Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Material</b>	Titanium

### MICRO STRAIGHT LOCKING PLATE 2.0MM

Product Code: B0.526

<b>Holes</b>	4 Hole, 5 Hole, 6 Hole, 7 Hole
<b>Locking Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Cortical Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Material</b>	Titanium



### MICRO RECON LOCKING PLATE 2.0MM

Product Code: **B0.527**

<b>Holes</b>	4 Hole To 12 Hole
<b>Locking Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Cortical Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Material</b>	Titanium



### MICRO T LOCKING PLATE 2.0MM

Product Code: **B0.528**

<b>Holes</b>	3 Hole, 5 Hole, 7 Hole
<b>Locking Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Cortical Screw 2.0mm Star Hex</b>	10mm to 16mm (Diff. 2.0mm)
<b>Material</b>	Titanium

## 2.4 MM LOCKING SCREW

Product Code: **B0.186**

<b>Size</b>	2.4mm
<b>Length</b>	06mm to 40mm (Diff. 2.0mm)
<b>Thread Diameter</b>	2.4mm
<b>Use Drills Diameter</b>	2.0mm
<b>Screw Driver</b>	Hex
<b>Material</b>	Stainless steel, Titanium



## 2.7 MM LOCKING SCREW

Product Code: **B0.204**

<b>Size</b>	2.7mm
<b>Length</b>	06mm to 40mm (Diff. 2.0mm)
<b>Thread Diameter</b>	2.7mm
<b>Use Drills Diameter</b>	2.0mm
<b>Screw Driver</b>	Hex
<b>Material</b>	Stainless Steel, Titanium



## 3.5 MM LOCKING SCREW

Product Code: **B0.168**

<b>Size</b>	3.5mm
<b>Length</b>	10mm to 60mm (Diff. 2.0mm)
<b>Thread Diameter</b>	3.5mm
<b>Use Drills Diameter</b>	2.7mm
<b>Screw Driver</b>	Hex
<b>Material</b>	Stainless Steel, Titanium



## 4.0 MM CANCELLOUS LOCKING SCREW

Product Code: **B0.170**

<b>Size</b>	4.0mm
<b>Length</b>	20mm to 85mm (Diff. 5.0mm)
<b>Thread Diameter</b>	4.0mm
<b>Use Drills Diameter</b>	3.2mm
<b>Screw Driver</b>	Hex
<b>Material</b>	Stainless steel, Titanium



## 5.0 MM LOCKING SCREW

Product Code: **B0.169**

<b>Size</b>	5.0mm
<b>Length</b>	16mm to 50mm (Diff. 2.0mm) 55mm to 95mm (Diff. 5.0mm)
<b>Thread Diameter</b>	5.0mm
<b>Use Drills Diameter</b>	4.0mm
<b>Screw Driver</b>	Hex
<b>Material</b>	Stainless Steel, Titanium

## 5.0 MM CANCELLOUS LOCKING SCREW

Product Code: **B0.9930**

<b>Size</b>	5.0mm
<b>Length</b>	30mm to 95mm (Diff. 5.0mm)
<b>Thread Diameter</b>	5.0mm
<b>Use Drills Diameter</b>	4.0mm
<b>Screw Driver</b>	Hex
<b>Material</b>	Stainless Steel, Titanium



## 6.5 MM LOCKING SCREW

Product Code: **B0.171**

<b>Size</b>	6.5mm
<b>Length</b>	30mm to 100mm (Diff. 5.0mm)
<b>Thread Diameter</b>	6.5mm
<b>Use Drills Diameter</b>	4.5mm
<b>Screw Driver</b>	Hex
<b>Material</b>	Stainless steel, Titanium



## 6.0MM LOCKING SCREW CANNULATED

Product Code: **B0.01.9930**

<b>Size</b>	Length
<b>6.0mm</b>	50mm to 120mm (Diff. 5.0mm)
<b>Thread Diameter</b>	6.0mm
<b>Drive</b>	Hex
<b>Material</b>	Stainless steel, Titanium

## 2.0 MM LOCKING SCREW (STAR)

Product Code: **B0.529**

<b>Size</b>	2.0mm
<b>Length</b>	06mm to 16mm (Diff. 2.0mm)
<b>Thread Diameter</b>	2.0mm
<b>Use Drills Diameter</b>	1.5mm
<b>Screw Driver</b>	Star
<b>Material</b>	Titanium



## WIRE MOUNT

Product Code: **B0.204**

<b>Size</b>	3.5mm,5.0mm,6.5mm
<b>Material</b>	Stainless Steel, Titanium





**Address:** No 64, Tingkat 1, Jalan Cempaka 1,  
Persiaran Cempaka, Bandar Aman Jaya,  
08000 Sungai Petani Kedah (Malaysia)

**Phone:** +6012 3433115 / +06 04 4205555/5556

**Email:** [mail@sahatechnology.com.my](mailto:mail@sahatechnology.com.my)

**Web:** [www.sahatechnology.com.my](http://www.sahatechnology.com.my)



TRAUMA  
IMPLANTS

 **SAHA**<sup>®</sup>  
Orthopedic Implants

Manufacturer of  
**Orthopedic** Implants &  
Surgical **Instruments**

## STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 MM

Product Code: **BOI.206**

<b>Holes</b>	5 Hole, 6 Hole, till 20 Hole
<b>Cortical Screw 3.5mm</b>	10mm to 60mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium



## CURVED PELVIC RECONSTRUCTION PLATE 3.5 MM

Product Code: **BOI.205**

<b>Holes</b>	6 Hole, 8 Hole, till 16 Hole
<b>Cortical Screw 3.5mm</b>	10mm to 60mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium

## PELVIC SPRING PLATE 3.5 MM

Product Code: **BOI.207**

<b>Holes</b>	1 Hole, 2 Hole, 3 Hole
<b>Cortical Screw 3.5mm</b>	10mm to 60mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium



**PELVIC RECONSTRUCTION J PLATE 3.5MM  
LEFT & RIGHT**

Product Code: **B0.208**

<b>Holes</b>	10Hole, 12Hole, 14Hole, 16Hole
<b>Locking Screw 3.5mm</b>	10mm to 60mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium



**PELVIC SYMPHYSEAL PLATE DC HOLE**

Product Code: **B0I.203**

<b>Holes</b>	4 Hole, 6 Hole, 8 Hole
<b>Cortical Screw 3.5mm</b>	10mm to 60mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium

## ACRO PLATE FOR 3.5 MM SCREW – LEFT & RIGHT

Product Code: **B0.9931**

<b>Holes</b>	4H,5H, 6H,7H
<b>Cortical Screw 3.5mm</b>	10mm to 60mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## LC DCP PLATE NARROW FOR 4.5 MM SCREW

Product Code: **B0I.134**

<b>Holes</b>	6H,7H,8H,9H,10H,11H,12H,13H,14
<b>Cortical Screw 4.5mm</b>	16mm to 70mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium

## DHS PLATE (DCP HOLES) SHORT BARREL

Product Code: **B0I.139**

<b>Holes</b>	2Hole To 18Hole
<b>Degree</b>	120° , 125° , 130° , 135° , 140° , 145°
<b>Lag Screw</b>	50mm to 115mm (Diff. 5.0mm)
<b>Cortical Screw 4.5mm</b>	16mm to 70mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## DHS PLATE(DOUBLE ANGLE)125D

Product Code: **BOI.01.140.125**

<b>Holes</b>	6Hole, 8Hole, 10Hole
<b>Lag Screw</b>	50mm to 115mm (Diff. 5.0mm)
<b>Cortical Screw 4.5mm</b>	16mm to 70mm (Diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium



## DHS PLATE (DCP HOLES) LONG BARREL

Product Code: **BOI.140**

<b>Holes</b>	2Hole To 18Hole
<b>Degree</b>	120°, 125°, 130°, 135°, 140°, 145°
<b>Lag Screw</b>	50mm to 115mm (Diff. 5.0mm)
<b>Cortical Screw 4.5mm</b>	16mm to 70mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium

## DHS PLATE PEDIATRICS

Product Code: **BO.931**

<b>Holes</b>	2H,3H,4H,5H,6H,7H,8H
<b>Degree</b>	125°, 130°, 135°
<b>Lag Screw</b>	50mm to 115mm (Diff. 5.0mm)
<b>Cortical Screw 4.5mm</b>	16mm to 70mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## 95° DCS PLATE (DCP HOLES)

Product Code: **B01.141**

<b>Holes</b>	6Hole To 20Hole
<b>Lag Screw</b>	50mm to 115mm (Diff. 5.0mm)
<b>Cortical Screw 4.5mm</b>	16mm to 70mm (Diff. 2.0mm)
<b>Cancellous Screw 6.5mm</b>	30mm to 110mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## TROCHANTERIC STABILIZATION PLATE

Product Code: **B0.9948**

<b>Holes</b>	4Hole
<b>Cortical Screw 4.5mm</b>	16mm to 70mm (Diff. 2.0mm)
<b>Cancellous Screw 6.5mm</b>	30mm to 110mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel

## EIGHT PLATE

Product Code: **B0.266**

<b>Size</b>	26mm, 29mm
<b>4.0mm Cannulated Screw</b>	12mm to 30mm (Diff. 2.0mm)
	30mm to 70mm (Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## EIGHT(HINGE) PLATE

Product Code: **BO.02.266**

<b>Size</b>	26mm, 29mm
<b>4.0mm Cannulated Screw</b>	12mm to 30mm (Diff. 2.0mm)
	30mm to 70mm (Diff. 5.0mm)
<b>Material</b>	Stainless steel, Titanium



## NARROW DYNAMIC COMPRESSION PLATE (DCP)

Product Code: **BOI.137**

<b>Holes</b>	6Hole To 14Hole
<b>Cortical Screw 4.5mm</b>	16mm to 70mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium

## OSTEOTOMI 'T' PLATE

Product Code: **BOI.187**

<b>Wedge Size</b>	7.5mm, 10mm, 12.5mm, 15mm
<b>Holes</b>	3Hole, 4Hole
<b>Cortical Screw 4.5mm</b>	16mm to 70mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## OSTEOTOMI PLATE (BUTTERFLY TYPE)

Product Code: **B0I.189**

<b>Wedge Size</b>	7.5mm,10mm,12.5mm,15mm
<b>Holes</b>	4Hole
<b>Cortical Screw 4.5mm</b>	16mm to 70mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## SMALL DYNAMIC COMPRESSION PLATE (DCP)

Product Code: **B0I.125**

<b>Holes</b>	3Hole to 16Hole
<b>Cortical Screw 3.5mm</b>	10mm to 40mm (Diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium

## 2.4MM CORTICAL SCREW

Product Code: **B0I.185**

<b>Size</b>	2.4mm
<b>Length</b>	06mm to 40mm (Diff. 2.0mm)
<b>Thread Diameter</b>	2.4mm
<b>Use Drills Diameter</b>	1.8mm
<b>Screw Driver</b>	2.4mm Hex
<b>Material</b>	Stainless steel, Titanium



## 2.7MM CORTICAL SCREW

Product Code: **B0I.101**

<b>Size</b>	2.7mm
<b>Length</b>	06mm to 40mm (Diff. 2.0mm)
<b>Thread Diameter</b>	2.7mm
<b>Use Drills Diameter</b>	2.2mm
<b>Screw Driver</b>	2.7mm Hex
<b>Material</b>	Stainless Steel, Titanium



## 3.5MM CORTICAL SCREW

Product Code: **B0I.103**

<b>Size</b>	3.5mm
<b>Length</b>	10mm to 60mm (Diff. 2.0mm)
<b>Thread Diameter</b>	3.5mm
<b>Use Drills Diameter</b>	2.5mm
<b>Screw Driver</b>	3.5mm Hex
<b>Material</b>	Stainless Steel, Titanium



## 4.0MM CANCELLOUS SCREW (SHORT THREAD)

Product Code: **B0I.104**

<b>Size</b>	4.0mm
<b>Length</b>	10mm to 30mm(Diff. 2.0mm) 35mm to 75mm(Diff. 5.0mm)
<b>Thread Diameter</b>	4.0mm
<b>Use Drills Diameter</b>	3.2mm
<b>Screw Driver</b>	3.5mm Hex
<b>Material</b>	Stainless steel, Titanium



## 4.0MM CANCELLOUS SCREW (FULL THREAD)

Product Code: **B0I.105**

<b>Size</b>	4.0mm
<b>Length</b>	10mm to 30mm(Diff. 2.0mm) 35mm to 75mm(Diff. 5.0mm)
<b>Thread Diameter</b>	4.0mm
<b>Use Drills Diameter</b>	3.2mm
<b>Screw Driver</b>	3.5mm Hex
<b>Material</b>	Stainless Steel, Titanium



## 4.0MM CANNULATED SCREW (SHORT THREAD)

Product Code: **B0I.106**

<b>Size</b>	4.0mm
<b>Length</b>	10mm to 30mm(Diff. 2.0mm) 35mm to 75mm(Diff. 5.0mm)
<b>Thread Diameter</b>	4.0mm
<b>Use Drills Diameter</b>	3.2mm
<b>Screw Driver</b>	3.5mm Hex
<b>Material</b>	Stainless Steel, Titanium



## 4.0MM CANNULATED SCREW (FULL THREAD)

Product Code: **B0I.107**


<b>Size</b>	4.0mm
<b>Length</b>	10mm to 30mm(Diff. 2.0mm) 35mm to 75mm(Diff. 5.0mm)
<b>Thread Diameter</b>	4.0mm
<b>Use Drills Diameter</b>	3.2mm
<b>Screw Driver</b>	3.5mm Hex
<b>Material</b>	Stainless steel, Titanium



## 4.5MM CORTICAL SCREW

Product Code: **B0I.109**

<b>Size</b>	4.5mm
<b>Length</b>	16mm to 70mm(Diff. 2.0mm)
<b>Thread Diameter</b>	4.5mm
<b>Use Drills Diameter</b>	3.2mm/3.5mm
<b>Screw Driver</b>	4.5mm Hex
<b>Material</b>	Stainless Steel, Titanium



## 4.5MM MALLEOLAR SCREW

Product Code: **B0I.110**

<b>Size</b>	4.5mm
<b>Length</b>	20mm to 100mm(Diff. 5.0mm)
<b>Thread Diameter</b>	4.5mm
<b>Use Drills Diameter</b>	3.2mm/3.5mm
<b>Screw Driver</b>	4.5mm Hex
<b>Material</b>	Stainless Steel, Titanium



## 6.5MM CANCELLOUS SCREW (16MM THREAD LENGTH)

Product Code: **B0I.111**

<b>Size</b>	6.5mm
<b>Length</b>	20mm to 110mm(Diff. 5.0mm)
<b>Thread Diameter</b>	6.5mm
<b>Use Drills Diameter</b>	4.0mm/4.5mm
<b>Screw Driver</b>	4.5mm Hex
<b>Material</b>	Stainless steel, Titanium



## 6.5MM CANCELLOUS SCREW (32MM THREAD LENGTH)

Product Code: **B0I.112**

<b>Size</b>	6.5mm
<b>Length</b>	20mm to 110mm(Diff. 5.0mm)
<b>Thread Diameter</b>	6.5mm
<b>Use Drills Diameter</b>	4.0mm/4.5mm
<b>Screw Driver</b>	4.5mm Hex
<b>Material</b>	Stainless Steel, Titanium



## 6.5MM CANCELLOUS SCREW (FULL THREAD)

Product Code: **B0I.113**

<b>Size</b>	6.5mm
<b>Length</b>	20mm to 110mm(Diff. 5.0mm)
<b>Thread Diameter</b>	6.5mm
<b>Use Drills Diameter</b>	4.0mm/4.5mm
<b>Screw Driver</b>	4.5mm Hex
<b>Material</b>	Stainless Steel, Titanium



## 6.5MM CANNULATED SCREW (16MM THREAD LENGTH)

Product Code: **B0I.114**

<b>Size</b>	6.5mm
<b>Length</b>	20mm to 110mm(Diff. 5.0mm)
<b>Thread Diameter</b>	6.5mm
<b>Use Drills Diameter</b>	4.0mm/4.5mm
<b>Screw Driver</b>	4.5mm Hex
<b>Material</b>	Stainless steel, Titanium



## 6.5MM CANNULATED SCREW (32MM THREAD LENGTH)

Product Code: **B0I.115**

<b>Size</b>	6.5mm
<b>Length</b>	20mm to 110mm(Diff. 5.0mm)
<b>Thread Diameter</b>	6.5mm
<b>Use Drills Diameter</b>	4.0mm/4.5mm
<b>Screw Driver</b>	4.5mm Hex
<b>Material</b>	Stainless Steel, Titanium



## 6.5MM CANNULATED SCREW (FULL THREAD)

Product Code: **B0I.117**

<b>Size</b>	6.5mm
<b>Length</b>	20mm to 110mm(Diff. 5.0mm)
<b>Thread Diameter</b>	6.5mm
<b>Use Drills Diameter</b>	4.0mm/4.5mm
<b>Screw Driver</b>	4.5mm Hex
<b>Material</b>	Stainless Steel, Titanium



## HERBERT SCREW

Product Code: **BOI.116**

<b>Size</b>	4.0mm
<b>Length</b>	10mm to 50mm(Diff. 2.0mm)
<b>Thread Diameter</b>	3.5mm/4.0mm & 4.0mm/3.0mm
<b>Use Drills Diameter</b>	2.7mm/3.2mm
<b>Screw Driver</b>	3.5mm Hex
<b>Material</b>	Stainless steel, Titanium



## HERBERT SCREW (HEADLESS)

Product Code: **BOI.119**

<b>Size</b>	Length
<b>3.0mm, 3.5mm, 4.0mm, 4.5mm, 5.5mm</b>	16mm to 80mm
<b>Material</b>	Titanium

## DHS LAG SCREW

Product Code: **BOI.118**

<b>Holes</b>	12.5mm
<b>Length</b>	50mm to 115mm(Diff. 5.0mm)
<b>Thread Diameter</b>	12.5mm
<b>Material</b>	Stainless Steel, Titanium



## DHS LAG SCREW (HEXAGONAL)

Product Code: **B0I.01.118**

<b>Size</b>	12.5mm
<b>Length</b>	50mm to 115mm(Diff. 5.0mm)
<b>Thread Diameter</b>	12.5mm
<b>Material</b>	Stainless Steel, Titanium



## DHS LAG SCREW PEDIATRICS

Product Code: **B0.263**

<b>Length</b>	40mm to 110mm(Diff. 5.0mm)
<b>Thread Diameter</b>	10mm
<b>Material</b>	Stainless steel, Titanium



**Address:** No 64, Tingkat 1, Jalan Cempaka 1,  
Persiaran Cempaka, Bandar Aman Jaya,  
08000 Sungai Petani Kedah (Malaysia)

**Phone:** +6012 3433115 / +06 04 4205555/5556

**Email:** [mail@sahatechnology.com.my](mailto:mail@sahatechnology.com.my)

**Web:** [www.sahatechnology.com.my](http://www.sahatechnology.com.my)



SPINE

 **SAHA**<sup>®</sup>  
Orthopedic Implants

Manufacturer of  
**Orthopedic** Implants &  
Surgical **Instruments**

## MONO AXIAL SCREW (DLX)

Product Code: **B0I.143**

<b>Diameter</b>	4.5mm,5.5mm,6.5mm,7.5mm
<b>Length</b>	25mm to 55mm(Diff. 5.0mm)
<b>Material</b>	Titanium



## POLY AXIAL SCREW (DLX)

Product Code: **B0I.144**

<b>Diameter</b>	4.5mm,5.5mm,6.5mm,7.5mm
<b>Length</b>	25mm to 55mm(Diff. 5.0mm)
<b>Material</b>	Titanium



## MONO AXIAL REDUCTION (DLX)

Product Code: **B0I.181**

<b>Diameter</b>	4.5mm,5.5mm,6.5mm,7.5mm
<b>Length</b>	25mm to 55mm(Diff. 5.0mm)
<b>Material</b>	Titanium

## POLY AXIAL REDUCTION (DLX)

Product Code: **B0I.182**

<b>Diameter</b>	4.5mm,5.5mm,6.5mm,7.5mm
<b>Length</b>	25mm to 55mm(Diff. 5.0mm)
<b>Material</b>	Titanium



## MINI MONO AXIAL SCREW

Product Code: **B0I.130**

<b>Diameter</b>	3.5mm,4.2mm
<b>Length</b>	12mm to 40mm(Diff. 2.0mm)
<b>Material</b>	Titanium

## MINI POLY AXIAL SCREW

Product Code: **B0I.131**

<b>Diameter</b>	3.5mm,4.2mm
<b>Length</b>	12mm to 40mm(Diff. 2.0mm)
<b>Material</b>	Titanium



## LAMINAR HOOK ANGLE BLADE

Product Code: B0.135.01

<b>Diameter</b>	Standard
<b>Length</b>	Standard
<b>Material</b>	Titanium



## LAMINAR HOOK NARROW BLADE

Product Code: B0.135 .02

<b>Diameter</b>	Standard
<b>Length</b>	Standard
<b>Material</b>	Titanium

## LEMINAR HOOK WIDE BLADE

Product Code: B0.135.03

<b>Diameter</b>	Standard
<b>Length</b>	Standard
<b>Material</b>	Titanium



## PEDICLE HOOK

Product Code: B0.135.04

<b>Diameter</b>	Standard
<b>Length</b>	Standard
<b>Material</b>	Titanium



## STAPLE

Product Code: B0.129

<b>Diameter</b>	Standard
<b>Length</b>	Standard
<b>Material</b>	Titanium

## TRANSVERSE CONNECTOR

Product Code: B0.138

<b>Diameter</b>	Standard
<b>Length</b>	Standard
<b>Material</b>	Titanium



## CONNECTING ROD

Product Code: **B0I.232**

<b>Diameter</b>	3.0mm, 3.5mm, 5.0mm, 5.5mm
<b>Length</b>	50mm to 120mm (Diff. 5.0mm)
	150mm to 500mm (Diff. 50.0mm)
<b>Material</b>	Titanium



## POSTERIOR OCCIPITAL PLATE (BUTTERFLY TYPE)


Product Code: **B0.132**

<b>Diameter</b>	3.5mm
<b>Length</b>	45mm,50mm,55mm,60mm
<b>Material</b>	Titanium

## POSTERIOR OCCIPITAL ROD

Product Code: **B0.123.3.5**

<b>Diameter</b>	3.5mm
<b>Length</b>	Standard
<b>Material</b>	Titanium



## POSTERIOR OCCIPITAL PLATE

Product Code: **B0.254**

<b>Diameter</b>	3.5mm
<b>Length</b>	200mm, 250mm, 300mm
<b>Material</b>	Titanium



## POSTERIOR OCCIPITAL SCREW

Product Code: **B0.133**

<b>Diameter</b>	4.0mm
<b>Length</b>	10mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Titanium

## ANTERIOR CERVICAL CAGE (SQUARE TYPE)

Product Code: **B0.144**

<b>Diameter</b>	4mm, To 10mm
<b>Length</b>	12mmx14mm, 12mmx16mm
<b>Material</b>	Titanium



## ANTERIOR CERVICAL PLATE - II SYS

Product Code: **B0.149**

<b>Diameter</b>	23mm to 74mm (Diff. 3.0mm)
<b>4mm Bone Screw</b>	12mm to 30mm
<b>Material</b>	Titanium



## CERVICLE CAGE OCTA

Product Code: **B0.513**

<b>Diameter</b>	3.5mm, 4.5mm, 5.5mm, 6.5mm, 7.5mm, 8.5mm
<b>Length</b>	12mm, 14mm, 16mm
<b>Material</b>	Titanium

## EXPANDABLE CAGE

Product Code: **B0.146**

<b>Diameter</b>	10mm, 12mm, 14mm, 16mm, 18mm
<b>Length</b>	25mm to 35mm 35mm to 55mm 55mm to 75mm
<b>Material</b>	Titanium



## ECO CAGE

Product Code: **B0.145**

<b>Diameter</b>	10mm, 12mm, 14mm, 16mm, 18mm, 20mm
<b>Length</b>	15mm to 100mm (Diff. 5.0mm)
<b>Material</b>	Titanium



## PLIF CAGE

Product Code: **B0.147**

<b>Diameter</b>	7mm, 8mm, 9mm, 10mm, 11mm
<b>Length</b>	20mm, 25mm, 30mm
<b>Material</b>	Titanium

## TLIF CAGE (BANANA TYPE)

Product Code: **B0.256**

<b>Diameter</b>	7mm, 8mm, 9mm, 10mm, 11mm
<b>Degree</b>	4°, 8°, 11°
<b>Material</b>	Titanium



## PLIF CAGE (BULLET TYPE)

Product Code: **B0.257**

<b>Diameter</b>	8mm, 10mm, 12mm, 14mm
<b>Length</b>	20mm, 25mm, 30mm
<b>Material</b>	Titanium



## ANTERIOR CERVICAL SCREW 4.0MM

Product Code: **B0.143**

<b>Diameter</b>	4.0mm
<b>Length</b>	10mm to 30mm (Diff. 2.0mm)
<b>Material</b>	Titanium

## OCCIPITAL LOOPS

Product Code: **B0.255**

<b>Diameter</b>	4mm
<b>Material</b>	Titanium






**Address:** No 64, Tingkat 1, Jalan Cempaka 1,  
Persiaran Cempaka, Bandar Aman Jaya,  
08000 Sungai Petani Kedah (Malaysia)

**Phone:** +6012 3433115 / +06 04 4205555/5556

**Email:** [mail@sahatechnology.com.my](mailto:mail@sahatechnology.com.my)

**Web:** [www.sahatechnology.com.my](http://www.sahatechnology.com.my)

**INTERLOCKING  
NAILINGS**



 **SAHA**<sup>®</sup>  
Orthopedic Implants

Manufacturer of  
**Orthopedic** Implants &  
Surgical **Instruments**

## SUPRA CONDYLAR NAIL

Product Code: **B0I.148**

<b>Nail Diameter</b>	9mm,10mm,11mm,12mm
<b>Nail Length</b>	15cm to 30cm(diff.5.0cm)
<b>6.5mm Bolt</b>	30mm to 100mm
<b>4.9mm Bolt</b>	22mm to 90mm
<b>Material</b>	Stainless steel, Titanium



## HUMERUS INTERLOCKING CANNULATED NAIL

Product Code: **B0I.149**

<b>Nail Diameter</b>	6mm,7mm,8mm
<b>Nail Length</b>	18cm to 30cm (diff. 2.0cm)
<b>Locking Bolt 3.1mm</b>	10mm to 50mm(diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## INTERLOCKING TIBIA NAIL

Product Code: **B0I.150**

<b>Nail Diameter</b>	8mm,9mm,10mm,11mm
<b>Nail Length</b>	28.0cm to 40.0cm(diff. 2.0cm)
<b>Locking Bolt 4.9mm</b>	20mm to 60mm(diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium

## INTERLOCKING FEMUR NAIL

Product Code: **B0I.151**

<b>Nail Diameter</b>	8mm,9mm,10mm,11mm,12mm
<b>Nail Length</b>	32.0cm to 46.0cm (Diff. 2.0cm)
<b>Locking Bolt 4.9mm</b>	20mm to 60mm(diff. 2.0mm)
<b>Material</b>	Stainless steel, Titanium



## TROCHANTERIC FEMUR NAIL (TFN)


Product Code: **B0.9935**

<b>Nail Diameter</b>	9mm,10mm,11mm,12mm
<b>Nail Length</b>	18cm
<b>Degree</b>	130° /135°
<b>Bolt 8.0mm</b>	40mm to 120mm(diff. 5.0mm)
<b>Bolt 6.4mm</b>	35mm to 120mm(diff. 5.0mm)
<b>Interlocking Bolt 4.9mm</b>	20mm to 60mm(diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium

## PROXIMAL FEMUR NAIL (PFN)

Product Code: **B0I.97**

<b>Nail Diameter</b>	9mm,10mm,11mm,12mm
<b>Nail Length</b>	25cm
<b>Degree</b>	130° /135°
<b>Bolt 8.0mm</b>	40mm to 120mm(diff. 5.0mm)
<b>Bolt 6.4mm</b>	35mm to 120mm(diff. 5.0mm)
<b>Interlocking Bolt 4.9mm</b>	20mm to 60mm(diff. 2.0mm)
<b>Material</b>	Stainless Steel, Titanium



## LONG PROXIMAL FEMUR NAIL - LEFT & RIGHT

Product Code: **B0I.97**

<b>Nail Diameter</b>	9mm,10mm,11mm,12mm
<b>Nail Length</b>	34 to 44cm (diff 2.0cm)
<b>Degree</b>	130° /135°
<b>Bolt 8.0mm</b>	40mm to 120mm(diff. 5.0mm)
<b>Bolt 6.4mm</b>	35mm to 120mm(diff. 5.0mm)
<b>Interlocking Bolt 4.9mm</b>	20mm to 60mm(diff. 2.0mm)
<b>Standard</b>	Left & Right
<b>Material</b>	Stainless steel, Titanium



## ELASTIC NAIL (FLEXIBLE NAIL)


Product Code: **B0I.9938**

<b>Nail Diameter</b>	2mm,2.5mm,3.0mm,3.5mm,4.0mm
<b>Nail Length</b>	44.0cm Standard
<b>Material</b>	Stainless Steel, Titanium

## ELASTIC NAIL (THREADED)

Product Code: **B0.01.9938**

<b>Nail Diameter</b>	2mm, 2.5mm, 3.0mm
<b>Nail Length</b>	18cm to 30cm (Diff. 2.0cm)
<b>Material</b>	Titanium



## SYRUS FEMUR NAIL - LEFT & RIGHT

Product Code: B01.159

<b>Nail Diameter</b>	9mm,10mm,11mm,12mm
<b>Nail Length</b>	34cm to 44cm (diff. 2.0cm)
<b>Cannulated 6.5mm Bolt</b>	30mm to 110m (diff. 5.0mm)
<b>Interlocking Bolt 4.9mm</b>	20mm to 60mm(diff. 2.0mm) 65mm to 80mm(diff. 5.0mm)
<b>Standard</b>	Left & Right
<b>Material</b>	Stainless Steel, Titanium



## PFNA-II PROXIMAL FEMUR NAIL ANTIROTATION (SHORT)

Product Code: B01.98

<b>Degree</b>	130°
<b>Nail Diameter</b>	9mm,10mm,11mm,12mm
<b>Nail Length</b>	17cm, 18cm,20cm,22cm,24cm
<b>Blade Size</b>	70mm to 115m (diff. 5.0mm)
<b>Interlocking Bolt 4.9mm</b>	20mm to 50mm(diff. 2.0mm)
<b>End Cap C.S</b>	5mm,10mm,15mm
<b>End Cap C.C</b>	5mm,10mm,15mm
<b>Material</b>	Stainless steel, Titanium



## PFNA-II PROXIMAL FEMUR NAIL ANTIROTATION (LONG)

Product Code: B01.98

<b>Degree</b>	130°
<b>Nail Diameter</b>	9mm,10mm,11mm,12mm
<b>Nail Length</b>	32cm to 44cm (diff. 2.0cm)
<b>Blade Size</b>	70mm to 115m (diff. 5.0mm)
<b>Interlocking Bolt 4.9mm</b>	20mm to 50mm(diff. 2.0mm)
<b>Standard</b>	Left & Right
<b>End Cap C.S</b>	5mm,10mm,15mm
<b>End Cap C.C</b>	5mm,10mm,15mm
<b>Material</b>	Stainless steel, Titanium



## SMART TIBIA INTERLOCKING NAIL

Product Code: **BOI.271**

<b>Nail Diameter</b>	9mm,10mm,11mm
<b>Nail Length</b>	28cm to 38cm (diff. 2.0cm)
<b>Interlocking Bolt 4.9mm</b>	20mm to 60mm(diff. 2.0mm)
	65mm to 80mm(diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## INTERLOCKING TIBIA NAIL(ANKLE)

Product Code: **BOI.9936**

<b>Nail Diameter</b>	9mm, 10mm, 11mm
<b>Nail Length</b>	15cm, 18cm, 21cm, 24cm, 27cm, 30cm, 33cm
<b>Material</b>	Titanium

## INTERLOCKING BOLT 3.1MM

Product Code: B0I.233

<b>Size</b>	3.1mm
<b>Length</b>	10mm to 50mm (Diff. 2.0mm)
<b>Thread Diameter</b>	3.1mm
<b>Use drills Diameter</b>	2.5mm
<b>Material</b>	Stainless steel, Titanium



## INTERLOCKING BOLT 4.9MM

Product Code: B0I.152

<b>Size</b>	4.9mm
<b>Length</b>	20mm to 60mm (Diff. 2.0mm) 65mm to 80mm (Diff. 5.0mm)
<b>Thread Diameter</b>	4.9mm
<b>Use drills Diameter</b>	4.0mm
<b>Material</b>	Stainless steel, Titanium

## PFN BOLT CANNULATED 6.4MM

Product Code: B0I.95

<b>Size</b>	6.4mm
<b>Length</b>	35mm to 120mm (Diff. 5.0mm)
<b>Thread Diameter</b>	6.4mm
<b>Material</b>	Stainless steel, Titanium



## PFN BOLT CANNULATED 8.0 MM

Product Code: **B0I.96**

<b>Size</b>	8.0mm
<b>Length</b>	35mm to 120mm (Diff. 5.0mm)
<b>Thread Diameter</b>	8mm
<b>Material</b>	Stainless steel, Titanium



## PFN BOLT CANNULATED 8.0MM (HEX)

Product Code: **B0I.01.196**

<b>Size</b>	8.0mm
<b>Length</b>	35mm to 120mm (Diff. 5.0mm)
<b>Thread Diameter</b>	8mm
<b>Material</b>	Stainless steel, Titanium

## PFNA II BLADE

Product Code: **B0I.89**

<b>Size</b>	70mm to 120mm
<b>Material</b>	Stainless steel, Titanium





**Address:** No 64, Tingkat 1, Jalan Cempaka 1,  
Persiaran Cempaka, Bandar Aman Jaya,  
08000 Sungai Petani Kedah (Malaysia)

**Phone:** +6012 3433115 / +06 04 4205555/5556

**Email:** [mail@sahatechnology.com.my](mailto:mail@sahatechnology.com.my)

**Web:** [www.sahatechnology.com.my](http://www.sahatechnology.com.my)



# HIP PROTHESIS

 **SAHA**<sup>®</sup>  
Orthopedic Implants

Manufacturer of  
**Orthopedic Implants &  
Surgical Instruments**

## THOMPSON HIP PROSTHESIS

Product Code: **B0I.154**

<b>Size</b>	37mm,39mm,41mm,43mm,45mm, 47mm,49mm,51mm,53mm
<b>Material</b>	Stainless steel



## AUSTIN MOORE HIP PROSTHESIS

Product Code: **B0I.155**

<b>Size</b>	37mm,39mm,41mm,43mm,45mm, 47mm,49mm,51mm,53mm
<b>Material</b>	Stainless steel

## BIPOLAR HIP PROSTHESIS

Product Code: **B0.9952**

<b>Size</b>	37mm,39mm,41mm,43mm,45mm, 47mm,49mm,51mm,53mm
<b>Material</b>	UHMWPE, Stainless steel



## MODULAR BIPOLAR CUP

Product Code: **B0.9950**

<b>Dia</b>	37mm, 39mm, 41mm, 43mm, 45mm,
	47mm, 49mm, 51mm, 53mm
<b>Size</b>	Small, Medium, Large
<b>Material</b>	UHMWPE Stainless steel



## MODULAR BIPOLAR STEM

Product Code: **B0.9951**

<b>Size</b>	Small, Medium, Large, XL, XXL
<b>Material</b>	Stainless steel

## SHOULDER PROSTHESIS

Product Code: **B0I.251**

<b>Size</b>	35, 37, 39, 43, 45, 47
<b>Material</b>	Stainless steel






**Address:** No 64, Tingkat 1, Jalan Cempaka 1,  
Persiaran Cempaka, Bandar Aman Jaya,  
08000 Sungai Petani Kedah (Malaysia)

**Phone:** +6012 3433115 / +06 04 4205555/5556

**Email:** [mail@sahatechnology.com.my](mailto:mail@sahatechnology.com.my)

**Web:** [www.sahatechnology.com.my](http://www.sahatechnology.com.my)

ARTHROSCOPY



 **SAHA**<sup>®</sup>  
Orthopedic Implants

Manufacturer of  
**Orthopedic** Implants &  
Surgical **Instruments**

## INTERFERENCE SCREW (ACL)

Product Code: **BOI.219**

<b>Diameter</b>	7.0mm, 8.0mm, 9.0mm, 10.0mm,11.0mm
<b>Screw Length</b>	20mm to 35mm(Diff. 5.0mm)
<b>Material</b>	Stainless Steel, Titanium



## ENDO BUTTON


Product Code: **BOI.218**

<b>Diameter</b>	Standard
<b>Length</b>	Standard
<b>Material</b>	Titanium

## ENDO BUTTON WITH LOOP

Product Code: **BOI.218**

<b>Size</b>	10mm, 15mm, 20mm, 25mm
<b>Material</b>	UHMWPE, Titanium



## SUTURE DISK

Product Code: **BOI.220**

<b>Diameter</b>	Standard
<b>4mm Bone Screw</b>	Standard
<b>Material</b>	Titanium



## SUTURE ANCHOR

Product Code: **BOI.217**

<b>Diameter</b>	3.0mm, 3.5mm, 4.0mm, 4.5mm, 5.0mm
<b>Length</b>	Standard
<b>Material</b>	Titanium

## SUTURE ANCHOR WITH LOOP

Product Code: **BOI.01.217**

<b>Diameter</b>	3.0mm, 5.0mm,
<b>Length</b>	Standard
<b>Material</b>	Titanium, UHMWPE



## LIGAMENT SPIKE WASHER(ROUND SHAPE)

Product Code: **B0I.254**

<b>Size</b>	Small, Medium, Large
<b>Material</b>	Titanium



## LIGAMENT SPIKE STAPLE(U SHAPE)

Product Code: **B0I.255**

<b>Size</b>	Small, Medium, Large
<b>Material</b>	Titanium



**Address:** No 64, Tingkat 1, Jalan Cempaka 1,  
Persiaran Cempaka, Bandar Aman Jaya,  
08000 Sungai Petani Kedah (Malaysia)

**Phone:** +6012 3433115 / +06 04 4205555/5556

**Email:** [mail@sahatechnology.com.my](mailto:mail@sahatechnology.com.my)

**Web:** [www.sahatechnology.com.my](http://www.sahatechnology.com.my)