


**INTERLOCKING
NAILINGS**



 **SAHA**[®]
Orthopedic Implants

Manufacturer of
Orthopedic Implants &
Surgical **Instruments**

SUPRA CONDYLAR NAIL

Product Code: **B0I.148**

Nail Diameter	9mm,10mm,11mm,12mm
Nail Length	15cm to 30cm(diff.5.0CM)
6.5mm Bolt	30mm to 100mm
4.9mm Bolt	22mm to 90mm
Material	Stainless steel, Titanium



HUMERUS INTERLOCKING CANNULATED NAIL

Product Code: **B0I.149**

Nail Diameter	6mm,7mm,8mm
Nail Length	18cm to 30cm (diff. 2.0cm)
Locking Bolt 3.1mm	10mm to 50mm(diff. 2.0mm)
Material	Stainless Steel, Titanium

INTERLOCKING TIBIA NAIL

Product Code: **B0I.150**

Nail Diameter	8mm,9mm,10mm,11mm
Nail Length	28.0cm to 40.0cm(diff. 2.0cm)
Locking Bolt 4.9mm	20mm to 60mm(diff. 2.0mm)
Material	Stainless Steel, Titanium



INTERLOCKING FEMUR NAIL

Product Code: **B0I.151**

Nail Diameter	8mm,9mm,10mm,11mm,12mm
Nail Length	32.0cm to 46.0cm (Diff. 2.0cm)
Locking Bolt 4.9mm	20mm to 60mm(diff. 2.0mm)
Material	Stainless steel, Titanium



TROCHANTERIC FEMUR NAIL (TFN)

Product Code: **B0.9935**

Nail Diameter	9mm,10mm,11mm,12mm
Nail Length	18cm
Degree	130° /135°
Bolt 8.0mm	40mm to 120mm(diff. 5.0mm)
Bolt 6.4mm	35mm to 120mm(diff. 5.0mm)
Interlocking Bolt 4.9mm	20mm to 60mm(diff. 2.0mm)
Material	Stainless Steel, Titanium

PROXIMAL FEMUR NAIL (PFN)

Product Code: **B0I.97**

Nail Diameter	9mm,10mm,11mm,12mm
Nail Length	25cm
Degree	130° /135°
Bolt 8.0mm	40mm to 120mm(diff. 5.0mm)
Bolt 6.4mm	35mm to 120mm(diff. 5.0mm)
Interlocking Bolt 4.9mm	20mm to 60mm(diff. 2.0mm)
Material	Stainless Steel, Titanium



LONG PROXIMAL FEMUR NAIL - LEFT & RIGHT

Product Code: **B01.97**

Nail Diameter	9mm,10mm,11mm,12mm
Nail Length	34 to 44cm (diff 2.0cm)
Degree	130° /135°
Bolt 8.0mm	40mm to 120mm(diff. 5.0mm)
Bolt 6.4mm	35mm to 120mm(diff. 5.0mm)
Interlocking Bolt 4.9mm	20mm to 60mm(diff. 2.0mm)
Standard	Left & Right
Material	Stainless steel, Titanium



ELASTIC NAIL (FLEXIBLE NAIL)

Product Code: **B01.9938**

Nail Diameter	2mm,2.5mm,3.0mm,3.5mm,4.0mm
Nail Length	44.0cm Standard
Material	Stainless Steel, Titanium

ELASTIC NAIL (THREADED)

Product Code: **B0.01.9938**


Nail Diameter	2mm, 2.5mm, 3.0mm
Nail Length	18cm to 30cm (Diff. 2.0cm)
Material	Titanium



SYRUS FEMUR NAIL - LEFT & RIGHT

Product Code: B01.159

Nail Diameter	9mm,10mm,11mm,12mm
Nail Length	34cm to 44cm (diff. 2.0cm)
Cannulated 6.5mm Bolt	30mm to 110m (diff. 5.0mm)
Interlocking Bolt 4.9mm	20mm to 60mm(diff. 2.0mm) 65mm to 80mm(diff. 5.0mm)
Standard	Left & Right
Material	Stainless Steel, Titanium



PFNA-II PROXIMAL FEMUR NAIL ANTIROTATION (SHORT)

Product Code: B01.98

Degree	130°
Nail Diameter	9mm,10mm,11mm,12mm
Nail Length	17cm, 18cm,20cm,22cm,24cm
Blade Size	70mm to 115m (diff. 5.0mm)
Interlocking Bolt 4.9mm	20mm to 50mm(diff. 2.0mm)
End Cap C.S	5mm,10mm,15mm
End Cap C.C	5mm,10mm,15mm
Material	Stainless steel, Titanium



PFNA-II PROXIMAL FEMUR NAIL ANTIROTATION (LONG)

Product Code: B01.98

Degree	130°
Nail Diameter	9mm,10mm,11mm,12mm
Nail Length	32cm to 44cm (diff. 2.0cm)
Blade Size	70mm to 115m (diff. 5.0mm)
Interlocking Bolt 4.9mm	20mm to 50mm(diff. 2.0mm)
Standard	Left & Right
End Cap C.S	5mm,10mm,15mm
End Cap C.C	5mm,10mm,15mm
Material	Stainless steel, Titanium



SMART TIBIA INTERLOCKING NAIL

Product Code: **BOI.271**

Nail Diameter	9mm,10mm,11mm
Nail Length	28cm to 38cm (diff. 2.0cm)
Interlocking Bolt 4.9mm	20mm to 60mm(diff. 2.0mm)
	65mm to 80mm(diff. 5.0mm)
Material	Stainless Steel, Titanium



INTERLOCKING TIBIA NAIL(ANKLE)

Product Code: **BOI.9936**

Nail Diameter	9mm, 10mm, 11mm
Nail Length	15cm, 18cm, 21cm, 24cm, 27cm, 30cm, 33cm
Material	Titanium

INTERLOCKING BOLT 3.1MM

Product Code: B0I.233

Size	3.1mm
Length	10mm to 50mm (Diff. 2.0mm)
Thread Diameter	3.1mm
Use drills Diameter	2.5mm
Material	Stainless steel, Titanium



INTERLOCKING BOLT 4.9MM

Product Code: B0I.152

Size	4.9mm
Length	20mm to 60mm (Diff. 2.0mm) 65mm to 80mm (Diff. 5.0mm)
Thread Diameter	4.9mm
Use drills Diameter	4.0mm
Material	Stainless steel, Titanium

PFN BOLT CANNULATED 6.4MM

Product Code: B0I.95

Size	6.4mm
Length	35mm to 120mm (Diff. 5.0mm)
Thread Diameter	6.4mm
Material	Stainless steel, Titanium



PFN BOLT CANNULATED 8.0 MM

Product Code: **B0I.96**

Size	8.0mm
Length	35mm to 120mm (Diff. 5.0mm)
Thread Diameter	8mm
Material	Stainless steel, Titanium



PFN BOLT CANNULATED 8.0MM (HEX)

Product Code: **B0I.01.196**

Size	8.0mm
Length	35mm to 120mm (Diff. 5.0mm)
Thread Diameter	8mm
Material	Stainless steel, Titanium

PFNA II BLADE

Product Code: **B0I.89**

Size	70mm to 120mm
Material	Stainless steel, Titanium





Address: No 64, Tingkat 1, Jalan Cempaka 1,
Persiaran Cempaka, Bandar Aman Jaya,
08000 Sungai Petani Kedah (Malaysia)

Phone: +6012 3433115 / +06 04 4205555/5556

Email: mail@sahatechnology.com.my

Web: www.sahatechnology.com.my