



SPINE

 **SAHA**[®]
Orthopedic Implants

Manufacturer of
**Orthopedic Implants &
Surgical Instruments**

MONO AXIAL SCREW (DLX)

Product Code: **B0I.143**

Diameter	4.5mm,5.5mm,6.5mm,7.5mm
Length	25mm to 55mm(Diff. 5.0mm)
Material	Titanium



POLY AXIAL SCREW (DLX)

Product Code: **B0I.144**

Diameter	4.5mm,5.5mm,6.5mm,7.5mm
Length	25mm to 55mm(Diff. 5.0mm)
Material	Titanium

MONO AXIAL REDUCTION (DLX)

Product Code: **B0I.181**

Diameter	4.5mm,5.5mm,6.5mm,7.5mm
Length	25mm to 55mm(Diff. 5.0mm)
Material	Titanium



POLY AXIAL REDUCTION (DLX)

Product Code: **B0I.182**

Diameter	4.5mm,5.5mm,6.5mm,7.5mm
Length	25mm to 55mm(Diff. 5.0mm)
Material	Titanium



MINI MONO AXIAL SCREW

Product Code: **B0I.130**

Diameter	3.5mm,4.2mm
Length	12mm to 40mm(Diff. 2.0mm)
Material	Titanium

MINI POLY AXIAL SCREW

Product Code: **B0I.131**

Diameter	3.5mm,4.2mm
Length	12mm to 40mm(Diff. 2.0mm)
Material	Titanium



LAMINAR HOOK ANGLE BLADE

Product Code: B0.135.01

Diameter	Standard
Length	Standard
Material	Titanium



LAMINAR HOOK NARROW BLADE

Product Code: B0.135 .02

Diameter	Standard
Length	Standard
Material	Titanium

LEMINAR HOOK WIDE BLADE

Product Code: B0.135.03

Diameter	Standard
Length	Standard
Material	Titanium



PEDICLE HOOK

Product Code: B0.135.04

Diameter	Standard
Length	Standard
Material	Titanium



STAPLE

Product Code: B0.129

Diameter	Standard
Length	Standard
Material	Titanium

TRANSVERSE CONNECTOR

Product Code: B0.138

Diameter	Standard
Length	Standard
Material	Titanium



CONNECTING ROD

Product Code: **B0I.232**

Diameter	3.0mm, 3.5mm, 5.0mm, 5.5mm
Length	50mm to 120mm (Diff. 5.0mm)
	150mm to 500mm (Diff. 50.0mm)
Material	Titanium



POSTERIOR OCCIPITAL PLATE (BUTTERFLY TYPE)


Product Code: **B0.132**

Diameter	3.5mm
Length	45mm,50mm,55mm,60mm
Material	Titanium

POSTERIOR OCCIPITAL ROD

Product Code: **B0.123.3.5**

Diameter	3.5mm
Length	Standard
Material	Titanium



POSTERIOR OCCIPITAL PLATE

Product Code: B0.254

Diameter	3.5mm
Length	200mm, 250mm, 300mm
Material	Titanium



POSTERIOR OCCIPITAL SCREW

Product Code: B0.133

Diameter	4.0mm
Length	10mm to 30mm (Diff. 2.0mm)
Material	Titanium

ANTERIOR CERVICAL CAGE (SQUARE TYPE)

Product Code: B0.144

Diameter	4mm, To 10mm
Length	12mmx14mm, 12mmx16mm
Material	Titanium



ANTERIOR CERVICAL PLATE - II SYS

Product Code: **B0.149**

Diameter	23mm to 74mm (Diff. 3.0mm)
4mm Bone Screw	12mm to 30mm
Material	Titanium



CERVICLE CAGE OCTA

Product Code: **B0.513**

Diameter	3.5mm, 4.5mm, 5.5mm, 6.5mm, 7.5mm, 8.5mm
Length	12mm, 14mm, 16mm
Material	Titanium

EXPANDABLE CAGE

Product Code: **B0.146**

Diameter	10mm, 12mm, 14mm, 16mm, 18mm
Length	25mm to 35mm 35mm to 55mm 55mm to 75mm
Material	Titanium



ECO CAGE

Product Code: **B0.145**

Diameter	10mm,12mm,14mm, 16mm, 18mm,20mm
Length	15mm to 100mm (Diff. 5.0mm)
Material	Titanium



PLIF CAGE

Product Code: **B0.147**

Diameter	7mm,8mm,9m, 10mm, 11mm
Length	20mm ,25mm,30mm
Material	Titanium

TLIF CAGE (BANANA TYPE)

Product Code: **B0.256**

Diameter	7mm,8mm,9m, 10mm, 11mm
Degree	4° , 8° , 11°
Material	Titanium



PLIF CAGE (BULLET TYPE)

Product Code: **B0.257**

Diameter	8mm, 10mm, 12mm, 14mm
Length	20mm, 25mm, 30mm
Material	Titanium



ANTERIOR CERVICAL SCREW 4.0MM

Product Code: **B0.143**

Diameter	4.0mm
Length	10mm to 30mm (Diff. 2.0mm)
Material	Titanium

OCCIPITAL LOOPS

Product Code: **B0.255**

Diameter	4mm
Material	Titanium





Address: No 64, Tingkat 1, Jalan Cempaka 1,
Persiaran Cempaka, Bandar Aman Jaya,
08000 Sungai Petani Kedah (Malaysia)

Phone: +6012 3433115 / +06 04 4205555/5556

Email: mail@sahatechnology.com.my

Web: www.sahatechnology.com.my