



TRAUMA
IMPLANTS

 **SAHA**[®]
Orthopedic Implants

Manufacturer of
Orthopedic Implants &
Surgical **Instruments**

STRAIGHT PELVIC RECONSTRUCTION PLATE 3.5 MM

Product Code: **BOI.206**

Holes	5 Hole, 6 Hole, till 20 Hole
Cortical Screw 3.5mm	10mm to 60mm (Diff. 2.0mm)
Material	Stainless steel, Titanium



CURVED PELVIC RECONSTRUCTION PLATE 3.5 MM

Product Code: **BOI.205**

Holes	6 Hole, 8 Hole, till 16 Hole
Cortical Screw 3.5mm	10mm to 60mm (Diff. 2.0mm)
Material	Stainless steel, Titanium

PELVIC SPRING PLATE 3.5 MM

Product Code: **BOI.207**

Holes	1 Hole, 2 Hole, 3 Hole
Cortical Screw 3.5mm	10mm to 60mm (Diff. 2.0mm)
Material	Stainless steel, Titanium



**PELVIC RECONSTRUCTION J PLATE 3.5MM
LEFT & RIGHT**

Product Code: **B0.208**

Holes	10Hole, 12Hole, 14Hole, 16Hole
Locking Screw 3.5mm	10mm to 60mm (Diff. 2.0mm)
Material	Stainless steel, Titanium



PELVIC SYMPHYSEAL PLATE DC HOLE

Product Code: **B0I.203**

Holes	4 Hole, 6 Hole, 8 Hole
Cortical Screw 3.5mm	10mm to 60mm (Diff. 2.0mm)
Material	Stainless steel, Titanium

ACRO PLATE FOR 3.5 MM SCREW – LEFT & RIGHT

Product Code: **B0.9931**

Holes	4H,5H, 6H,7H
Cortical Screw 3.5mm	10mm to 60mm (Diff. 2.0mm)
Material	Stainless Steel, Titanium



LC DCP PLATE NARROW FOR 4.5 MM SCREW

Product Code: **B0I.134**

Holes	6H,7H,8H,9H,10H,11H,12H,13H,14
Cortical Screw 4.5mm	16mm to 70mm (Diff. 2.0mm)
Material	Stainless Steel, Titanium

DHS PLATE (DCP HOLES) SHORT BARREL

Product Code: **B0I.139**

Holes	2Hole To 18Hole
Degree	120° , 125° , 130° , 135° , 140° , 145°
Lag Screw	50mm to 115mm (Diff. 5.0mm)
Cortical Screw 4.5mm	16mm to 70mm (Diff. 2.0mm)
Material	Stainless Steel, Titanium



DHS PLATE(DOUBLE ANGLE)125D

Product Code: **BOI.01.140.125**

Holes	6Hole, 8Hole, 10Hole
Lag Screw	50mm to 115mm (Diff. 5.0mm)
Cortical Screw 4.5mm	16mm to 70mm (Diff. 2.0mm)
Material	Stainless steel, Titanium



DHS PLATE (DCP HOLES) LONG BARREL

Product Code: **BOI.140**

Holes	2Hole To 18Hole
Degree	120°, 125°, 130°, 135°, 140°, 145°
Lag Screw	50mm to 115mm (Diff. 5.0mm)
Cortical Screw 4.5mm	16mm to 70mm (Diff. 2.0mm)
Material	Stainless Steel, Titanium



DHS PLATE PEDIATRICS

Product Code: **BO.931**

Holes	2H,3H,4H,5H,6H,7H,8H
Degree	125°, 130°, 135°
Lag Screw	50mm to 115mm (Diff. 5.0mm)
Cortical Screw 4.5mm	16mm to 70mm (Diff. 2.0mm)
Material	Stainless Steel, Titanium



95° DCS PLATE (DCP HOLES)

Product Code: **B01.141**

Holes	6Hole To 20Hole
Lag Screw	50mm to 115mm (Diff. 5.0mm)
Cortical Screw 4.5mm	16mm to 70mm (Diff. 2.0mm)
Cancellous Screw 6.5mm	30mm to 110mm (Diff. 5.0mm)
Material	Stainless Steel, Titanium



TROCHANTERIC STABILIZATION PLATE

Product Code: **B0.9948**

Holes	4Hole
Cortical Screw 4.5mm	16mm to 70mm (Diff. 2.0mm)
Cancellous Screw 6.5mm	30mm to 110mm (Diff. 5.0mm)
Material	Stainless Steel

EIGHT PLATE

Product Code: **B0.266**

Size	26mm, 29mm
4.0mm Cannulated Screw	12mm to 30mm (Diff. 2.0mm)
	30mm to 70mm (Diff. 5.0mm)
Material	Stainless Steel, Titanium



EIGHT(HINGE) PLATE

Product Code: **BO.02.266**

Size	26mm, 29mm
4.0mm Cannulated Screw	12mm to 30mm (Diff. 2.0mm)
	30mm to 70mm (Diff. 5.0mm)
Material	Stainless steel, Titanium



NARROW DYNAMIC COMPRESSION PLATE (DCP)

Product Code: **BOI.137**

Holes	6Hole To 14Hole
Cortical Screw 4.5mm	16mm to 70mm (Diff. 2.0mm)
Material	Stainless Steel, Titanium

OSTEOTOMI 'T' PLATE

Product Code: **BOI.187**

Wedge Size	7.5mm, 10mm, 12.5mm, 15mm
Holes	3Hole, 4Hole
Cortical Screw 4.5mm	16mm to 70mm (Diff. 2.0mm)
Material	Stainless Steel, Titanium



OSTEOTOMI PLATE (BUTTERFLY TYPE)

Product Code: **B0I.189**

Wedge Size	7.5mm,10mm,12.5mm,15mm
Holes	4Hole
Cortical Screw 4.5mm	16mm to 70mm (Diff. 2.0mm)
Material	Stainless Steel, Titanium



SMALL DYNAMIC COMPRESSION PLATE (DCP)

Product Code: **B0I.125**

Holes	3Hole to 16Hole
Cortical Screw 3.5mm	10mm to 40mm (Diff. 2.0mm)
Material	Stainless Steel, Titanium

2.4MM CORTICAL SCREW

Product Code: **B0I.185**

Size	2.4mm
Length	06mm to 40mm (Diff. 2.0mm)
Thread Diameter	2.4mm
Use Drills Diameter	1.8mm
Screw Driver	2.4mm Hex
Material	Stainless steel, Titanium



2.7MM CORTICAL SCREW

Product Code: **B0I.101**

Size	2.7mm
Length	06mm to 40mm (Diff. 2.0mm)
Thread Diameter	2.7mm
Use Drills Diameter	2.2mm
Screw Driver	2.7mm Hex
Material	Stainless Steel, Titanium



3.5MM CORTICAL SCREW

Product Code: **B0I.103**

Size	3.5mm
Length	10mm to 60mm (Diff. 2.0mm)
Thread Diameter	3.5mm
Use Drills Diameter	2.5mm
Screw Driver	3.5mm Hex
Material	Stainless Steel, Titanium



4.0MM CANCELLOUS SCREW (SHORT THREAD)

Product Code: **B0I.104**

Size	4.0mm
Length	10mm to 30mm(Diff. 2.0mm) 35mm to 75mm(Diff. 5.0mm)
Thread Diameter	4.0mm
Use Drills Diameter	3.2mm
Screw Driver	3.5mm Hex
Material	Stainless steel, Titanium



4.0MM CANCELLOUS SCREW (FULL THREAD)

Product Code: **B0I.105**

Size	4.0mm
Length	10mm to 30mm(Diff. 2.0mm) 35mm to 75mm(Diff. 5.0mm)
Thread Diameter	4.0mm
Use Drills Diameter	3.2mm
Screw Driver	3.5mm Hex
Material	Stainless Steel, Titanium



4.0MM CANNULATED SCREW (SHORT THREAD)

Product Code: **B0I.106**

Size	4.0mm
Length	10mm to 30mm(Diff. 2.0mm) 35mm to 75mm(Diff. 5.0mm)
Thread Diameter	4.0mm
Use Drills Diameter	3.2mm
Screw Driver	3.5mm Hex
Material	Stainless Steel, Titanium



4.0MM CANNULATED SCREW (FULL THREAD)

Product Code: **B0I.107**


Size	4.0mm
Length	10mm to 30mm(Diff. 2.0mm) 35mm to 75mm(Diff. 5.0mm)
Thread Diameter	4.0mm
Use Drills Diameter	3.2mm
Screw Driver	3.5mm Hex
Material	Stainless steel, Titanium



4.5MM CORTICAL SCREW

Product Code: **B0I.109**

Size	4.5mm
Length	16mm to 70mm(Diff. 2.0mm)
Thread Diameter	4.5mm
Use Drills Diameter	3.2mm/3.5mm
Screw Driver	4.5mm Hex
Material	Stainless Steel, Titanium



4.5MM MALLEOLAR SCREW

Product Code: **B0I.110**

Size	4.5mm
Length	20mm to 100mm(Diff. 5.0mm)
Thread Diameter	4.5mm
Use Drills Diameter	3.2mm/3.5mm
Screw Driver	4.5mm Hex
Material	Stainless Steel, Titanium



6.5MM CANCELLOUS SCREW (16MM THREAD LENGTH)

Product Code: **B0I.111**

Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)
Thread Diameter	6.5mm
Use Drills Diameter	4.0mm/4.5mm
Screw Driver	4.5mm Hex
Material	Stainless steel, Titanium



6.5MM CANCELLOUS SCREW (32MM THREAD LENGTH)

Product Code: **B0I.112**

Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)
Thread Diameter	6.5mm
Use Drills Diameter	4.0mm/4.5mm
Screw Driver	4.5mm Hex
Material	Stainless Steel, Titanium



6.5MM CANCELLOUS SCREW (FULL THREAD)

Product Code: **B0I.113**

Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)
Thread Diameter	6.5mm
Use Drills Diameter	4.0mm/4.5mm
Screw Driver	4.5mm Hex
Material	Stainless Steel, Titanium



6.5MM CANNULATED SCREW (16MM THREAD LENGTH)

Product Code: **B0I.114**

Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)
Thread Diameter	6.5mm
Use Drills Diameter	4.0mm/4.5mm
Screw Driver	4.5mm Hex
Material	Stainless steel, Titanium



6.5MM CANNULATED SCREW (32MM THREAD LENGTH)

Product Code: **B0I.115**

Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)
Thread Diameter	6.5mm
Use Drills Diameter	4.0mm/4.5mm
Screw Driver	4.5mm Hex
Material	Stainless Steel, Titanium



6.5MM CANNULATED SCREW (FULL THREAD)

Product Code: **B0I.117**

Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)
Thread Diameter	6.5mm
Use Drills Diameter	4.0mm/4.5mm
Screw Driver	4.5mm Hex
Material	Stainless Steel, Titanium



HERBERT SCREW

Product Code: **BOI.116**

Size	4.0mm
Length	10mm to 50mm(Diff. 2.0mm)
Thread Diameter	3.5mm/4.0mm & 4.0mm/3.0mm
Use Drills Diameter	2.7mm/3.2mm
Screw Driver	3.5mm Hex
Material	Stainless steel, Titanium



HERBERT SCREW (HEADLESS)

Product Code: **BOI.119**

Size	Length
3.0mm, 3.5mm, 4.0mm, 4.5mm, 5.5mm	16mm to 80mm
Material	Titanium

DHS LAG SCREW

Product Code: **BOI.118**

Holes	12.5mm
Length	50mm to 115mm(Diff. 5.0mm)
Thread Diameter	12.5mm
Material	Stainless Steel, Titanium



DHS LAG SCREW (HEXAGONAL)

Product Code: **B0I.01.118**

Size	12.5mm
Length	50mm to 115mm(Diff. 5.0mm)
Thread Diameter	12.5mm
Material	Stainless Steel, Titanium



DHS LAG SCREW PEDIATRICS

Product Code: **B0.263**

Length	40mm to 110mm(Diff. 5.0mm)
Thread Diameter	10mm
Material	Stainless steel, Titanium



Address: No 64, Tingkat 1, Jalan Cempaka 1,
Persiaran Cempaka, Bandar Aman Jaya,
08000 Sungai Petani Kedah (Malaysia)

Phone: +6012 3433115 / +06 04 4205555/5556

Email: mail@sahatechnology.com.my

Web: www.sahatechnology.com.my